

SEQUENCE LISTING

<110> Roche Diagnostics GmbH  
F. Hoffmann-La Roche AG

<120> Improved detection of bacterial (mollicutes) contamination

<130> 26237

<150> EP09009965.6  
<151> 2009-08-01

<160> 4

<170> PatentIn version 3.5

<210> 1  
<211> 21  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> universal primer (forward) for the detection of Mycoplasma and  
related species

<400> 1  
ggcgaatggg tgagtaacac g 21

<210> 2  
<211> 22  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> universal primer (reverse) for the detection of Mycoplasma and  
related species

<400> 2  
cggataacgc ttgcgaccta tg 22

<210> 3  
<211> 21  
<212> DNA  
<213> Chinese hamster ovary (CHO) cells

<220>  
<221> misc\_feature  
<223> forward primer specific for a target sequence in the  
Glyceraldehyde 3-phosphate dehydrogenase gene

<400> 3  
caaaggcaca gtcaaggctg a 21

<210> 4  
<211> 21  
<212> DNA  
<213> Chinese hamster ovary (CHO) cells

<220>  
<221> misc\_feature  
<223> reverse primer specific for a target sequence in the  
Glyceraldehyde 3-phosphate dehydrogenase gene

<400> 4  
tgggtgaagac gccagtagat t